



Ag/AgCl

6.0733.100

c(KCl) = 3 mol/L/
14/15

6.0733.100

Qt.	Order no.	
1 PCS	6.2008.040	Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Order no.**6.2106.020****Strand / 1 m / 2 x B**

Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B

**6.2106.060****Strand / 2 m / 2 x B**

Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B

**6.2308.020****Electrolyte 3 mol/L KCl (250 mL)**

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)

