



Hamilton PRP-X100 - 250/4.0


6.1005.010

Hamilton PRP-X100 - 250/4.0 IC

Order no.

6.02821.000 **Holder for the PRP-X100 Guard cartridge**

Holder for the Hamilton PRP-X 100 Guard cartridge for connection to the analytical separation column.



6.1005.020

PRP-X100 guard column cartridge/4.0

For protecting Hamilton PRP-X100 columns. The cartridge efficiently removes contamination in the form of particles due, for example, to the growth of bacteria or algae.

