

Metrosep Dual 4 - 100/4.6

6.1016.030

Metrosep Dual 45 mL/Min

Metrosep Dual 4 - 100/4.693 g/L0.5 μ g/LEPA 314

6.1016.030

Qt.	Order no.
2 PCS	6.2744.060 Threaded stopper
	For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Order no.

6.1016.500 Metrosep Dual 4 Guard Column Kit/4.6

Even if the Dual-4 columns based on monolithic silica gel are very durable, the use of the Dual 4 guard column is advised in order to further increase the safety of the analytical separation column. The Dual 4 guard column is a PEEK cartridge which is also filled with monolithic silica gel. This cartridge is easy to replace and is screwed directly onto the analytical column in an aluminum holder. The proven "On Column Guard System" is simple to use and also offers the advantage of minimal dead volume.



6.1016.510 Metrosep Dual 4 Guard Column Cartridge/4.6

Set with three guard column cartridges

