



# Metrosep Carb 2 - 100/4.0

6.1090.410

ICMetrosep Carb 2 - 100/4.0/ pH = 014 Metrosep Carb 2100mm

## 6.1090.410

Qt.	Order no.	
2 PCS	<b>6.2744.060</b>	<b>Threaded stopper</b>
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.		



**Order no.****6.1015.200      Metrosep BO3 Trap 1 - 100/4.0**

Trap column for the removal of borate contaminants from the eluent. The Metrosep BO<sub>3</sub> Trap 1 - 100/4.0 is used mainly in carbohydrate analysis with hydroxide eluents.

**6.1015.300      Metrosep CO3 Trap 1 - 100/4.0**

Trap column for the removal of carbonate traces in hydroxide eluents.

**6.1090.500      Metrosep Carb 2 Guard/4.0**

The Metrosep Carb 2 Guard/4.0 effectively removes contamination, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.

**6.1090.510      Metrosep Carb 2 S-Guard/4.0**

The Metrosep Carb 2 S-Guard/4.0 effectively removes contaminations, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.

