




Multi-Mode Electrode pro

6.1246.120

Mercury electrode for voltammetry. Can be operated as DME, SMDE or HMDE.

Order no.

6.1226.030 **Mercury drop capillary (non-silanized) for MME**

Non-silanized mercury drop capillary (10 pieces)



6.1226.050 Mercury drop capillary (silanized) for MME

Silanized mercury drop capillary (10 pieces)



6.1247.020 Needles for MME

Needles for Multi Mode Electrode (3 pieces).



6.1247.110 Screw cap for the mercury refill opening for MME pro



6.1247.120 Retaining screw for capillary for MME pro



6.1247.140

Needle holder for MME pro



6.1247.150

Silicone seal, green 13.9 mm for MME pro



6.1247.210

Assembly ring for needle holder for MME



6.1247.220

Wrench for MME



6.2112.130

Adapter plug contact 5.5 mm / screw contact, gold-plated

For VA electrodes



6.2615.130

Electrode stand for MME pro

Electrode stand for filling and storing the Multi-Mode Electrode pro.

