



Vision Air 2.0

6.6072.201

Vision Air -

Vision Air

Vision Air :

-
-
- SQL

Vision Air (66072201) :

-
-

VisionPLS_Toolbox Unscrambler®

Order no.**6.6072.300 Pre-calibration for the lumber industry**

Pre-calibration for the determination of the kappa number and density as well as breaking, tensile and buckling strength in wood pulp by means of NIR spectroscopy.

**6.6072.301 Pre-calibration for the analysis of polyols**

Pre-calibration for determination of the hydroxyl number in polyols by means of NIR spectroscopy.

**6.6072.302 Pre-calibration for crude palm oil analysis**

Determination of iodine number, moisture, and free fatty acids in crude palm oil.



6.6072.303

Pre-calibration for crude palm kernel oil analysis

Determination of iodine number, moisture, and free fatty acids in crude palm kernel oil.



6.6072.304

Pre-calibration for RBD palm oil analysis

Determination of iodine number, moisture, and free fatty acids in RBD palm oil



6.6072.305

Pre-calibration for diesel analysis

Determination of cetane index und cetane number, as well as density, cold filter plugging point, flash point and sulfur in diesel fuel.



6.6072.306

Pre-calibration for gasoline analysis


Determination of octane number, LCI, benzene, aromatic compounds, olefins, oxygen and density in gasoline.



6.6072.307

Pre-calibration for kerosine analysis

Determination of cetane index, aromatic compounds, density, flash point, hydrogen, viscosity and other parameters in kerosine.



6.6072.312

Pre-calibration for isocyanate

Pre-calibration for determination of the NCO number in isocyanate using NIR spectroscopy.

