


Solid Autosampler CIC (TEI)

2.0136.0620

The Solid Autosampler CIC (TEI) is used for automated analysis of solid samples using Combustion IC. It can be used in combination with the Boat Introduction Modules (6.07311.010 or 6.07311.020). The sample rack for solids no. 1 (6.07312.010) and a cover for the sample rack are already included.

Phạm vi giao hàng 2.0136.0620

Qt.	Order no.	Description
1 PCS	1.0136.0620	Solid Autosampler CIC (TEI)



1 PCS


6.07312.010

Solid sample tray no. 1

Solid sample tray no. 1 for the Solid Autosampler CIC (2.0136.0620) for 20 positions.



Phụ kiện tùy chọn

Order no.	Description	
6.07312.020	Solid sample tray no. 2	
	Solid sample tray no. 2 for the Solid Autosampler CIC (2.0136.0620) for 20 positions.	
6.07312.030	Solid sample tray no. 3	
	Solid sample tray no. 3 for the Solid Autosampler CIC (2.0136.0620) for 20 positions.	