

## Combined polymer membrane electrode, NO3

6.00510.120

Combined nitrate-selective electrode with polymer membrane. This ISE is suitable for:

- ion measurements of  $NO_3^-$  (7\*10<sup>-6</sup> to 1 mol/L) in aqueous solutions
- titrations in aqueous solutions

Thanks to a robust/break-proof plastic shaft made of propylene and an impact protection for the polymer membrane, this sensor is mechanically very resistant.

The reference electrolyte used is c(KCI) = 3 mol/L.

Phạm vi giao hàng 6.00510.120

Qt. Order no. Description



Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



## Phụ kiện tùy chọn

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F	
	For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL)	_
	Electrolyte solution $c(KCI) = 3 \text{ mol/L}$ , (for Ag/AgCl reference systems)	
		The state of the s
6.2308.050	Electrolyte 3 mol/L KCl (50 mL)	
	Electrolyte solution $c(KCI) = 3 \text{ mol/L (for Ag/AgCl reference systems)}$	
		The state of the s