



Swing arm for Polytron and 786 Swing Head, left swinging

6.1462.250



Swing arm for Polytron and 786 Swing Head, left swinging, with one buret tip for the addition of solvent or titrant and three spray nozzles for cleaning. The swing arm allows the usage of the Polytron for sample preparation on multi row racks.

Phạm vi giao hàng 6.1462.250

Qt.	Order no.	Description
1 PCS	6.1808.000	Connection bushing M6



Phụ kiện tùy chọn

Order no.	Description	
6.1829.050	Replacement tubings for retracted tubings in titration heads and robotic arms Three replacement tubings with tapering tip.	
6.9012.000	Dispersing aggregate for Polytron, 125 mm Standard dispersing aggregate for Polytron 1300 D	
6.9012.010	Dispersing aggregate for Polytron, 157 mm Special dispersing aggregate (Metrohm version) with protruding knives	