

IC equipment: Dosino regeneration

6.5330.190

Accessory set for the assembly of a Dosino for automatic regeneration of the Metrohm Suppressor Module (MSM).

Phạm vi giao hàng 6.5330.190

Qt.	Order no.	Description
1 PCS	6.1580.120	807 Dosing Unit 2 mL without accessories

807 Dosing Unit with integrated data chip with 2 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread.



1 PCS

6.1602.160

Bottle attachment GL 45 for eluents and reagents

For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.




1 PCS

6.1618.020

Thread adapter S 40 to GL 45

Adapter for bottles with S40 thread (Merck...) to GL45.



1 PCS

6.1619.000

Adsorber tube for 807 Dosing Unit



1 PCS

6.1624.000

Adapter SGJ 14 for 6.1619.XXX Adsorber tube

Adapter SGJ 14. Part of 6.1619.010 for Dilution Sample Processor.




1 PCS 6.1805.120 FEP tubing / M6 / 100 cm
With light and kink protection



1 PCS 6.1829.020 FEP aspiration tube M6, 0.5 m
For Dosing Unit



1 PCS 6.2057.210 Holder for Dosino to IC instruments
Mounting part for Dosinos at the IC instrument housing.



1 PCS 6.2744.080 M6 thread / UNF 10/32 adapter
With M6 outer thread. For IC instruments.



1 PCS

6.2821.120

Inline filter 2 µm

Inline filter 2 µm

