

# 833

2.833.0030

## Scope of delivery 2.833.0030

Qt.	Order no.	Description	
1 PCS	6.1803.020	PTFE capillary 0.97 mm i.d. / 5 m	
	For all IC instrum	ents	



Connecting IC Interface - IC Pump Unit, IC Suppressor Module, IC Dialysis Unit or IC Liquid Handling Unit



2 PCS

6.2744.010 Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

6.2744.030 Coupling olive/UNF 10/32, 4x

Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump



1 PCS

6.2744.040 2 x UNF 10/32 coupling

For connecting 1/16 in. capillaries. For IC instruments



Inline filter 2  $\mu m$ 





### **Optional accessories**

Order no.	Description		
6.1803.030	PTFE capillary 0.5 mm i.d. / 3 m		
	Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit		

#### 6.1826.060 Pump tubing (orange/yellow)

For Compact IC, Personal IC, Basic IC



#### 6.2755.000 Tubing cartridge

For the peristaltic pumps of IC instruments. Contact pressure adjustable.



#### 6.2821.130 Spare filter for inline filter

Spare filter for inline filter.





#### 6.2835.000 SPM Rotor A

Rotor for inline sample preparation with cation exchange



### 6.2835.010 Connecting piece to MSM and SPM module

Universal connecting piece for MSM in Professional IC Vario and Compact Flex instruments and for MSM-HC and SPM module. With inlets and outlets.

