



# Sample Processor 843 Pump Station– rinse /aspirate

2.843.0150

843 Pump Station Sample Processor /rinse/aspirate 814 815 Sample Processor

## Scope of delivery 2.843.0150

Qt.	Order no.	Description
1 PCS	1.843.0120	<b>843 Pump Station (peristaltic)</b>  The 843 Pump Station (peristaltic) has two built-in peristaltic pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.



1 PCS

**6.1543.170**

**Aspiration tip / M8 thread**

Aspiration tip for the Sample Processors.



---

2 PCS

**6.1621.000**

**Container 10 L**

As rinsing or waste container in automated systems.



---

2 PCS

**6.1621.100**

**Cap for 10-L container**

Cap for 10-L container. For easy and safe transportation of filled containers.



---

4 PCS

**6.1801.120**

**PVC tubing / 4 mm / 6 mm / 2 m**



3 PCS

**6.1805.060**

**FEP tubing / M6 / 60 cm**

With light and kink protection



---

3 PCS

**6.1805.420**

**FEP tubing / M6 / 48 cm**

With light and kink protection



---

1 PCS

**6.1805.510**

**PTFE tubing / M8 / 60 cm**

With protection against kinks.



---

1 PCS

**6.1808.170**

**Distributor for Sample Processors for rinsing  
and aspirating**

Distributor for installation on Sample Processors for rinsing and aspirating when using externally connected pumps.



2 PCS

**6.1812.000**

**PTFE tubing 4 mm / 6 mm / 4 m**



---

4 PCS

**6.1820.050**

**M6 / Olive adapter**

For 6.1826.100 pump tubing, to peristaltic pumps.



---

1 PCS

**6.1826.100**

**Pump tubing 6.4 mm / 1 m**

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



---

2 PCS

**6.1828.000**

**PVDF connection nipple**

For 6.1621.000 canister.



1 PCS

**6.2141.300**

**Remote cable 843 Pump Station to Sample Processor**

Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors



---

3 PCS


**6.2740.020**

**Spray nozzle**

For the fine-spraying of the rinsing solution.



# Optional accessories

Order no.	Description	
6.2740.030	<p><b>Rinsing nozzle</b></p> <p>The rinsing solution is applied as a fine jet for the optimal removal of layers on the electrodes and titration accessories. For 6.5613.000 ECO titration equipment.</p>	

---