



NIRS XDS RapidLiquid

2.921.0410

NIRS XDS RapidLiquid Analyzer NIRS XDS RapidLiquid Analyzer

Scope of delivery 2.921.0410




Qt.	Order no.	Description
-----	-----------	-------------

Rapid, precise analyses of liquids and suspensions of all types

The NIRS XDS RapidLiquid Analyzer enables rapid, precise analyses of liquid formulations and substances. Precise measurement results at the push of a button make the NIRS XDS RapidLiquid Analyzer an equally reliable and simple solution for quality monitoring in laboratories and processes. The samples are transferred to quartz cuvettes designed for multiple use or disposable glass vials; a tempered sample compartment ensures reproducible analysis conditions and thus accurate measurement results.



Optional accessories

Order no.	Description	
6.7402.070	Disposable vials, 2 mm diameter, transmission 200 lockable disposable glass vials (borosilicate) with a diameter of 2 mm for analyses of liquids in transmission mode. Compatible with: <ul style="list-style-type: none">• Holder OMNIS NIR, vial, 2 mm (6.07401.050)• DS2500 holder for flow cells (6.7492.000)	
6.7402.100	Disposable vials, 10 mm diameter, transmission 195 lockable disposable glass vials (borosilicate) with a diameter of 10 mm for analyses of liquids in transmission. Compatible with: <ul style="list-style-type: none">• Holder OMNIS NIR, vial, 10 mm (6.07401.080)• DS2500 holder for 4 mm disposable vials (6.7492.030)	
6.7402.240	Disposable vial, 8 mm, transmission, qty. 100 100 disposable glass vials (borosilicate) with an optical path length of 8 mm for analyses of liquids in transmission. The disposable vials are supplied with the associated stoppers (number of pieces = 100). Compatible with: <ul style="list-style-type: none">• Holder OMNIS NIR, vial, 8 mm (6.07401.070)• DS2500 holder for 8 mm disposable vials (6.7492.020)	

6.7402.260 Stopper for disposable vial, 4 mm, qty. 100

100 stoppers for 4 mm transmission disposable vials (6.7402.250).



6.7403.000 NIRS XDS spacer for disposable 8 mm vials (6.7402.000)

Spacer for disposable glass vials with 8 mm diameter.



6.7403.020 XDS spacer for disposable 2 mm vials (6.7402.070)

Spacer for the disposable glass vials with 2 mm diameter.



6.7403.030 XDS spacer for disposable 10 mm vials (6.7402.100)

Spacer for the disposable glass vials with 10 mm diameter.

