



dSolvotrode

6.00203.300

/ OMNIS pH

dTodes OMNIS Titratoren

Optional accessories

Order no.	Description
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m

Cable for connecting a dTode (plug-in head Q) to the digital measuring module from OMNIS.



6.02104.310

Electrode cable plug-in head Q / plug P, 1.5 m

Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.



6.2320.000

TEABr 0.4 mol/L in ethylene glycol (250 mL)

Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol), $c(\text{TEABr}) = 0.4 \text{ mol/L}$

