

**5** •

6.1414.030

-

## Scope of delivery 6.1414.030

| Qt.   | Order no.         | Description                        |  |
|-------|-------------------|------------------------------------|--|
| 1 PCS | 6.1454.130        | Form seal made of FPM for 61414xxx |  |
|       | for KF and VA ins | struments                          |  |



## **Optional accessories**

| Order no.  | Description  |
|------------|--|
| 6.1244.040 | Set of seals   |
|            | Set of seals for 6.1414.030. The set is comprised of one seal and various O-rings. |

## 6.2730.010 Screw nipple

Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations.



## 6.2730.020 Septum stopper

With E.301.0041 O-ring. For 6.1448.010 Septum and 6.1414.030 Titration vessel lid. Used for KF applications.

