



# 159 x 2 mL 3 x 300 mL

6.2041.480

858 Professional Sample Processor 159 3 2 mL (Eppendorf) 6.1608.0806.1627.000

(PVC)

## Optional accessories

Order no.	Description
6.1608.080	Bottle / 300 mL

Rinsing vessel for IC Sample Processors.



**6.1627.000**

**PE bottle, 125 mL for IC**

This bottle is used with fully automatic dilution.



---

**6.1627.100**

**Screw cap with septum, trace analysis**

The screw cap with septum is ideal for trace analysis due to the negligible contamination caused by the septum. Perfect for maintaining standards over a prolonged period. Screw cap with septum for PE bottle, 125 mL (6.1627.000) und PE bottle, 50 mL (6.1608.100).



---

**6.2072.000**

**Bottle adapter for 6.1627.000 to 6.2041.440**

The bottle adapter permits placement of the 125 mL PE bottle on the 148 x 11 mL + 3 x 300 mL (6.2041.440) sample rack.

