




MBA.S

MBA.S MultiBAMBAAutolab / MBA MBA 5 MBA.S FRA32M.S

Optional accessories

Order no.	Description
RRDE.AUPT.S	Rotating ring disk electrode, Gold disk Platinum ring
	The RRDE electrode tip consist of a 5 mm disc of gold with a concentric platinum ring at a distance of 375µm, leading to a theoretical collection efficiency of 24.9%. 

RRDE.GCPT.S Rotating ring disk electrode, Glassy Carbon disk Platinum ring

The RRDE electrode tip consist of a 5 mm disc of Glassy Carbon with a concentric platinum ring at a distance of 375 μ m, leading to a theoretical collection efficiency of 24.9%.

