



OMNIS Rod Stirrer Titration

2.1006.0020

Stabrührer zum Anschluss an ein OMNIS Titrator (Kabellänge 1.1 m).

Lieferumfang 2.1006.0020

Qt.	Order no.	Beschreibung
1 PCS	1.1006.0020	Rod Stirrer Titration
		Rod stirrer for connection to one OMNIS Titrator (cable length 1.1 m).



Optionales Zubehör

Order no.	Beschreibung	
6.01900.010	Stirring propeller 30 mm ETFE Stirring propeller 30 mm ETFE with adjustable height for use with Rod Stirrer "Titrator" or Rod Stirrer "Sample Robot"	
6.01900.030	Stirring propeller 20 mm ETFE Stirring propeller, 20 mm, with adjustable height for use with Rod Stirrer "Titrator" or Rod Stirrer "Sample Robot"	