



Nadeldrucker TX-900MH

2.143.0200

Kompakt-Nadeldrucker mit USB und RS232 Schnittstelle für

- 900 Touch Control (+)
- 915 KF Ti-Touch (+)
- 916 Ti-Touch (+)
- 917 Coulometer (+)
- 91X Meter (+)
- Eco Dosimat / Titrator (+)
- 877 / 848 Titrino plus (zusätzlich 6.2151.100 benötigt)
- 865 / 876 Dosimat plus (zusätzlich 6.2151.100 benötigt)
- 862 Compact Titrosampler (zusätzlich 6.2151.100 benötigt)
- 870 KF Titrino plus (zusätzlich 6.2151.100 benötigt)
- 899 Coulometer (zusätzlich 6.2151.100 benötigt)

(+) Kabel im Lieferumfang enthalten

Lieferumfang 2.143.0200

Qt.	Order no.	Beschreibung
-----	-----------	--------------

1 PCS	1.143.0200	TX-900MH impact printer
-------	------------	-------------------------

Compact impact printer with USB and RS232 interface for

- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877/848 Titrino plus
- 865/876 Dosimat plus
- 91X meter (cable 6.2151.140)
- Eco Dosimat
- Eco Titrator
- Eco Coulometer
- 862 Compact Titrosampler
- 870 KF Titrino plus
- 899 Coulometer



1 PCS	6.2134.040	RS-232 DB9 - IBM-PC(DB9) connecting cable
-------	------------	---

For the Metrohm RS-232 interface - PC connection. For Titrinos, Titrino plus and 756, 831 KF Coulometer, RS-232/USB Box



1 PCS	6.2151.020	Cable USB A - USB B 1.8 m
-------	------------	---------------------------

USB connecting cable



Y cable for connection of a USB printer to the pH/Conductometers 912 / 913 / 914. This cable allows a printer and the power supply unit to be connected to the measuring instrument at the same time.



Optionales Zubehör

Order no.	Beschreibung	
6.2151.100	Adapter USB MINI (OTG) - USB A For connecting USB instruments.	
6.9001.070	Paper roll for TX-900MH Paper roll for TX-900MH impact printer	
6.9001.080	Continuous adhesive labels for TX-900MH Self-adhesive paper roll for TX-900MH impact printer	
6.9002.070	Color ribbon for TX-900MH impact printer	