


IC-VA-Detektor 791

2.791.0010

D

Lieferumfang 2.791.0010

Qt.	Order no.	Beschreibung
1 PCS	1.791.0010	791 IC VA Detector
	D	



1 PCS	6.2813.020	Dummy Cell for VA Detector
		Electrical simulation of a voltammetry cell, used for simulation and control purposes. For VA Detectors



Optionales Zubehör

Order no.	Beschreibung	
6.1204.130	Silver Electrode Tip Metal electrode tip, silver (Ag). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.	
6.1204.140	Gold Electrode Tip Metal electrode tip, gold (Au). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.	