



800 Dosino

2.800.0010

Der 800 Dosino ist ein Antrieb mit Schreib-/Lesehardware für intelligente Dosiereinheiten. Mit fest montiertem Kabel (Länge 150 cm).

Lieferumfang 2.800.0010

Qt.	Order no.	Beschreibung
1 PCS	1.800.0030	800 Dosino (cable length 1.5 m)

Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 1.5 m).



1 PCS

6.1546.030

Piston tongs

For pistons for 807 Dosing Units or cylinder units OMNIS.



Optionales Zubehör

Order no.	Beschreibung	
6.1562.130	Transfer Tubing with holder 2 x M6 / 10 mL to Sample Processors Transfer tubing for dilution and preconcentration using a Dosino and Sample Processor.	
6.2151.010	Cable MSB plug / socket Extension cable for MSB connections and for Touch Control - Titrande, Mini-DIN socket - mini-DIN plug	
6.3032.120	807 Dosing Unit 2 mL 807 Dosing Unit with integrated data chip with 2 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.	
6.3032.150	807 Dosing Unit 5 mL 807 Dosing Unit with integrated data chip with 5 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.	

6.3032.210

807 Dosing Unit 10 mL

807 Dosing Unit with integrated data chip with 10 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



6.3032.220

807 Dosing Unit 20 mL

807 Dosing Unit with integrated data chip with 20 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



6.3032.250

807 Dosing Unit 50 mL

807 Dosing Unit with integrated data chip with 50 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.

