

# pH Lab Irda 120V US

2.827.0015

## Lieferumfang 2.827.0015

| Qt.                                                                                   | Order no.  | Beschreibung                                                               |
|---------------------------------------------------------------------------------------|------------|----------------------------------------------------------------------------|
| 1 PCS                                                                                 | 6.2013.010 | Clamping ring<br><br>For support rods with a diameter of 10 mm.            |
|    |            |                                                                            |
| 1 PCS                                                                                 | 6.2016.050 | Support rod / 300 mm                                                       |
|   |            |                                                                            |
| 1 PCS                                                                                 | 6.2021.020 | Electrode holder<br><br>Electrode holder for 4 electrodes and 2 buret tips |
|  |            |                                                                            |

1 PCS

6.2161.020

Plug-in power supply unit US 120V/6V DC

Power supply for 827 pH lab



---

1 PCS

6.2621.070

5 mm hex key for IC Sample Processors

