



827 pH lab; 230 V, EU mit Primatrode

2.827.0114

Labor-pH-Meter mit IrDa-Schnittstelle und Primatrode 6.0228.010

Lieferumfang 2.827.0114

| Qt. | Order no. | Beschreibung |
|-------|------------|---------------|
| 1 PCS | 6.2013.010 | Clamping ring |

For support rods with a diameter of 10 mm.



1 PCS

6.2016.050

Support rod / 300 mm



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2161.010

Plug-in power supply unit EU 230V /6V DC

Power supply for 827 pH lab





1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



Optionales Zubehör

| Order no. | Beschreibung | |
|------------|--|---|
| 6.2307.230 | Buffer solutions pH 4, 7 and 9 Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL |  |
| 6.2325.000 | pHit kit Maintenance kit for pH electrodes The kit contains: <ul style="list-style-type: none">• 50 mL cleaning solution• 50 mL 3M KCl solution• 50 mL storage solution• 2 Storage vessels• Instructions for use |  |