




900 Touch Control

2.900.0010

Bedienungseinheit für die Titrandos, USB Sample Processoren, 856 Conductivity Module, 867 pH Module und 846 Dosing Interface. Berührungssensitives, hochauflösendes Farbdisplay, einfache und intuitive Bedienung dank Favoriten für den direkten Methodenzugriff. Mit integrierter Ethernet-Schnittstelle für direkte Anbindung ans Internet und USB-Schnittstelle für den Anschluss von USB-Druckern oder eines USB-Memorysticks.

Dialogsprachen: Deutsch, Englisch, Chinesisch, Französisch, Spanisch, Portugiesisch, Russisch, Koreanisch, Polnisch und Italienisch.

Lieferumfang 2.900.0010


Qt.	Order no.	Beschreibung
-----	-----------	--------------

1 PCS

1.900.0010

900 Touch Control

Operating unit for the Titrando, USB Sample Processors, 856 Conductivity Module, 867 pH Module and 846 Dosing Interface. Touch-sensitive, high-resolution color display, simple and intuitive operation, thanks to Favorites for direct method access. With integrated Ethernet interface for direct connection to the Internet and USB interface for connecting USB printers or a USB memory stick.



Dialog languages: German, English, Chinese, French, Spanish, Portuguese, Russian, Korean, Polish and Italian.

1 PCS

6.2723.300

Protective cover for 900 Touch Control

Protective cover for 900 Touch Control



1 PCS

6.6064.010

USB Memory Stick for 900 Touch Control

