



909 UV Digester (230 V)

2.909.0014

Aufschlussgerät für die UV-Fotolyse von Wasserproben mit niedriger bis mittlerer organischer Belastung. Zur Probenvorbereitung in der Spurenelementbestimmung mittels Voltammetrie, Ionenchromatographie und Spektroskopie (AAS, ICP). Integriertes Gerät mit Bedienungseinheit und Nassteil. Mit Luftkühlung und automatischer Steuerung von Aufschlusstemperatur und -zeit. Für 12 Proben mit je maximal 12 mL Probevolumen. Gerät für 220 ... 240 V und 50 ... 60 Hz.

Lieferumfang 2.909.0014

Qt.	Order no.	Beschreibung
1 PCS	1.909.0014	909 UV Digester (230 V)

Digestion instrument for UV photolysis of water samples with low to medium organic load. For sample preparation in trace element determination by means of voltammetry, ion chromatography and spectroscopy (AAS, ICP). Integrated instrument with operating unit and wet end. With air cooling and automatic control of digestion temperature and time. For 12 samples with a maximum of 12 mL sample volume each. Instrument for 220 - 240 V and 50 - 60 Hz.




1 PCS

6.1110.010

Pt1000 temperature sensor for 909 UV Digester

Pt1000 temperature sensor (class B) made of glass with fixed cable (0.5 m) for connecting to a 909 UV Digester.



12 PCS

6.1446.100

PTFE stopper for 12 mL UV quartz sample vessel

For 6.2414.000 quartz sample vessel



1 PCS

6.1446.190

Thermometer stopper for Pt1000 temperature sensor

For 6.1110.010 Pt1000 temperature sensor and 6.2414.000 12 mL quartz sample vessel.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



1 PCS

6.2041.240

Sample holder for 12 quartz sample vessels 12 mL

For 909 UV Digester.



1 PCS

6.2225.000

UV warning plate

For 909 UV Digester



12 PCS

6.2414.000

Quartz sample vessel 12 mL

For UV Digester



1 PCS

6.2621.100

3 mm hex key for IC Sample Processors



1 PCS

6.2745.100

UV protective shield

For 909 UV Digester.



1 PCS

6.2804.090

UV mercury vapor lamp

Medium pressure mercury vapor lamp for 909 UV Digester.



Optionales Zubehör

Order no.	Beschreibung	
6.2042.050	Spacer for Pt1000 temperature sensor	
