



919 IC Autosampler plus for CVS

2.919.0120

Der 919 IC Autosampler plus for CVS eignet sich für kleine Probenserien vollautomatischer Bestimmungen von organischen Additiven in galvanischen Bädern wie Glanzbildner («Brightener») und Grundeinebner («Suppressoren») mittels «Cyclic Voltammetric Stripping» (CVS) zusammen mit Professional-VA/CSVs-Geräten und viva. Das System wird direkt von der viva-Software gesteuert und erlaubt den automatischen Methodenwechsel zwischen verschiedenen Bestimmungsarten. Das mitgelieferte Probenrack fasst 14 Proben von maximal 50 mL.

Der 919 IC Autosampler plus for CVS wird inklusive Zubehör, Probenrack und Anschlusskabel für den Aufbau eines kompletten Arbeitsplatzes ausgeliefert, jedoch ohne Spülausrüstung.

Lieferumfang 2.919.0120

Qt.	Order no.	Beschreibung
-----	-----------	--------------

1 PCS

1.919.0020

919 IC Autosampler plus

The 919 IC Autosampler plus fulfills the requirements of laboratories with medium sample numbers. It enables automation of the full range of Metrohm IC instruments.



1 PCS

6.1458.080

Transfer head for 919 IC Autosampler plus for VA and CVS applications.



1 PCS

6.1803.020

PTFE capillary 0.97 mm i.d. / 5 m

For all IC instruments



1 PCS

6.1815.010

Spiral band 0.5 m

For holding together different cables or tubing.




2 PCS

6.1826.020

PVC pump tubing (blue/blue)

For online IC instruments and automation in voltammetry.



1 PCS

6.1831.020

PEEK capillary 0.75 mm i.d., 22 cm

For all IC components.



1 PCS

6.1835.050

PEEK sample tube

For sample changers for transferring larger-sized sample volumes, e.g. for brightener determination.



1 PCS

6.2041.490

Sample rack, 28 x 50 mL, (PVC)

Sample rack for 919 Autosampler plus, 858 Professional Sample Processor suitable for 28 sample vessels. The following sample vessels can be used: 6.2747.010.

Plastic: Polyvinyl chloride (PVC)



1 PCS

6.2151.000

Cable USB A – mini-DIN 8-pin

Controller cable



1 PCS

6.2621.080

Capillary cutter

Cutting tool for simple cutting of plastic capillaries, 1/16" and 1/8", including 5 replacement blades.



1 PCS

6.2621.090

1/2 in. wrench

For 1/2 in. screws. For IC Sample Processors



1 PCS

6.2621.100

3 mm hex key for IC Sample Processors



1 PCS 6.2621.130 Hexagon key 2 mm
2 mm.



1 PCS 6.2739.000 Wrench
For tightening connectors



2 PCS 6.2744.010 Pressure screw 5x
With UNF 10/32 connection. For the connection of PEEK capillaries



2 PCS 6.2744.060 Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



1 PCS

6.2744.080

M6 thread / UNF 10/32 adapter

With M6 outer thread. For IC instruments.




2 PCS

6.2744.160

Nozzle/UNF 10/32 coupling with security device


Security device for connecting pump tubing to a nozzle.



1 PCS

6.2744.290

T connector 3 × UNF 10/32 with holder



4 PCS

6.2747.010

Sample vessel 50 mL

With blue screw cap and mL graduation. Set of 25 pieces.



1 PCS

6.2752.110

Holder for T connector to 919 IC Autosampler plus



1 PCS

6.2833.030

Needle holder 1/8 in. for 863 Compact Autosampler

Needle holder 1/8 in. for 863 Compact Autosampler

