



Ecotrode Gel mit Pt1000

6.00221.600

Wartungsfreie pH-Elektrode (Gelelektrolyt) mit integriertem Temperatursfühler (Pt1000) für Routine-pH-Messungen in gleichartigen Proben. Die Elektrode wird in c(KCl) = sat. (6.2308.000) gelagert und ist nicht geeignet für ionenarme Lösungen.

Der Alterungsindikator zeigt einen nötigen Austausch der Elektrode frühzeitig an.

Lieferumfang 6.00221.600


Qt.	Order no.	Beschreibung
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Optionales Zubehör

Order no.	Beschreibung	
6.2104.600	Electrode cable for plug in head U/plug F, 2x2 mm B, 1m For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2104.610	Electrode cable for plug in head U/plug F, 2x2 mm B, 2m For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2308.000	Electrolyte KCl (sat.), 250 mL Electrolyte KCl (sat.), 250 mL	