



Einstichelektrode mit Pt1000

6.00226.600

Wartungsfreie kombinierte pH-Elektrode (Gelelektrolyt) für Einstichmessungen aller Art (z.B. bei Käse, Fleisch, Teig) mit integriertem Pt1000-Temperaturfühler. Die Elektrode wird in gesättigter Kaliumchloridlösung $c(\text{KCl}) = \text{sat.}$ (6.2308.000) gelagert und ist nicht geeignet für ionenarme Lösungen. Der Alterungsindikator zeigt einen nötigen Austausch der Elektrode frühzeitig an.

Lieferumfang 6.00226.600


Qt.	Order no.	Beschreibung
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Optionales Zubehör

Order no.	Beschreibung	
6.2104.600	Electrode cable for plug in head U/plug F, 2x2 mm B, 1m For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2104.610	Electrode cable for plug in head U/plug F, 2x2 mm B, 2m For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2308.000	Electrolyte KCl (sat.), 250 mL Electrolyte KCl (sat.), 250 mL	