



iAg-Titrode

6.00470.300

Intelligente, kombinierte Silberringelektrode mit einer pH-Glasmembran als Referenzelektrode und integriertem Speicherchip für Sensordaten.

Diese wartungsfreie Elektrode eignet sich für Fällungstitrationsen bei konstant bleibendem pH-Wert (Titriermittel Silbernitrat) von z.B.:



- Chlorid, Bromid, Iodid
- Sulfiden
- Schwefelwasserstoff
- Mercaptane
- Cyanide

Diese Elektrode wird in dest. Wasser aufbewahrt.

Je nach Anwendung empfiehlt sich die Verwendung einer Ag-Titrode mit Ag_2S -Überzug, die entsprechend bestellt werden kann.

iTodes können an Titrando, Ti-Touch oder den 913/914 Metern verwendet werden.

Lieferumfang 6.00470.300

Qt.	Order no.	Beschreibung
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.
		
1 PCS	6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.
		

Optionales Zubehör

Order no.	Beschreibung
2.854.0010	<div>854 iConnect</div> <div>The 854 iConnect – Electrode cable and measuring amplifier for intelligent "iTrodes" electrodes (cable length 150 cm).</div>

