



## Leitfähigkeitsmesszelle $c = 0.1 \text{ cm}^{-1}$ mit Pt1000 (Fixkabel, 2.0 m)

6.00922.040

Leitfähigkeitsmesszelle aus Edelstahl mit Zellkonstante  $c = 0.1 \text{ cm}^{-1}$  (Richtwert) mit integriertem Pt1000-Temperaturfühler und Fixkabel (2.0 m) zur Anbindung an das OMNIS Messmodul Leitfähigkeit.

Dieser Sensor eignet sich für automatisierte Messungen von tiefen Leitfähigkeiten ( $0.01 \dots 1000 \mu\text{S}/\text{cm}$ ) wie z.B. in deion. Wasser oder für Messungen nach USP <645>.

Lieferumfang 6.00922.040

| Qt.   | Order no.  | Beschreibung                            |
|-------|------------|-----------------------------------------|
| 1 PCS | 6.1236.050 | Sleeve with SGJ 14/12 mm                |
|       |            | Sleeve with SGJ 14/12 mm, polyethylene. |



# Optionales Zubehör

| Order no. | Beschreibung |
|-----------|--------------|
|-----------|--------------|

6.02101.020      Measuring module conductivity

Measurement channel for one OMNIS Titrator or Titration Module for the connection of conductivity measuring cells.



6.2324.010      Conductivity standard 100 µS/cm 250 mL

Conductivity standard for calibration of conductivity measuring cells with cell constant = 0.1/cm.



6.2324.110      Conductivity standard 100 µS/cm, 5 x 30 mL

Conductivity standard for calibration of conductivity measuring cells with cell constant = 0.1/cm.

