



Leitfähigkeitsmesszelle $c = 0.5 \text{ cm}^{-1}$ mit Pt1000 (Fixkabel, 2.0 m)

6.00924.080

4-Leiter-Leitfähigkeitsmesszelle mit Zellkonstante $c = 0.5 \text{ cm}^{-1}$ (Richtwert), mit integriertem Pt1000-Temperaturfühler und Fixkabel (2.0 m) zur Anbindung an das OMNIS Messmodul Leitfähigkeit.

Dank einem robustem/bruchsicherem Kunststoffschacht aus PEEK ist dieser Sensor mechanisch sehr belastbar und eignet sich für automatisierte Messungen von mittleren Leitfähigkeiten ($0.01 \mu\text{S}/\text{cm}$ bis $200 \text{ mS}/\text{cm}$) wie z.B. in

- Trinkwasser
- Oberflächenwasser
- Abwasser

Lieferumfang 6.00924.080

Qt.	Order no.	Beschreibung
-----	-----------	--------------

1 PCS


6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



Optionales Zubehör

Order no.	Beschreibung
6.02101.020	Measuring module conductivity Measurement channel for one OMNIS Titrator or Titration Module for the connection of conductivity measuring cells. 
6.2301.060	KCl ion standard 250 mL KCl ion standard ($c(\text{ion}) = 0.1000 \pm 0.0005 \text{ mol/L}$), can also be used as standard solution for conductometry 12.87 mS/cm (25°C). 
6.2324.010	Conductivity standard $100 \mu\text{S/cm}$ 250 mL Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$. 
6.2324.110	Conductivity standard $100 \mu\text{S/cm}$, 5 x 30 mL Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$. 