



# Metrosep A Supp 18 - 250/4.0

6.01033.430

Die Selektivität der Metrosep A Supp 18 Säulen ist spezifisch auf die Arbeit mit Hydroxid-Eluenten ausgelegt. Die Metrosep A Supp 18 - 250/4.0 ist die längere A Supp 18 Säulenversion. Durch die kleine Partikelgrösse von 3.5 µm werden hohe Trennleistungen erreicht. Die symmetrischen Peaks erlauben die einfache Bestimmung von Anionen, auch in tiefen µg/L Konzentrationsbereichen. Die Metrosep A Supp 18 - 250/4.0 ist besonders für die Trennung der Oxohalogenide Chlorit, Bromat und Chlorat neben den Standardanionen bei Raumtemperatur geeignet.

## Optionales Zubehör

Order no.	Beschreibung
-----------	--------------

The Metrosep A Supp 18 Guard/4.0 reliably protects the Metrosep A Supp 18 analytical separation column against contamination from sample or eluent. It contains the same separation material as the Metrosep A Supp 18 IC Columns, is also made of PEEK, and is screwed directly onto the respective separating column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, without influencing chromatographic separating efficiency.

---

