



Metrosep C 6 - 100/2.0

6.01051.210

Die 100-mm-Version der microbore Metrosep-C-6-Säule ist für die Bestimmung der Standardkationen beispielsweise in Trinkwasser bestimmt. Es werden kurze Retentionszeiten erzielt, dies mit einer verhältnismässig guten Natrium/Ammonium-Trennung. Die hohe Kapazität des C-6-Materials erlaubt höhere Probenvolumen.

Die Säule ist für den Einsatz in der IC-MS-Kopplung geeignet.

Lieferumfang 6.01051.210

Qt.	Order no.	Beschreibung
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optionales Zubehör

Order no.	Beschreibung
6.01051.600	Metrosep C 6 Guard/2.0

The Metrosep C 6 Guard/2.0 contains the C 6 column material and is used to protect against particles and contamination. This considerably lengthens the service life of the analytical separation column. The Metrosep C 6 Guard/2.0 works according to the "On Column Guard System" and is screwed directly onto the respective separation column with virtually no dead volume.

