



Metrosep C Supp 2 - 100/4.0

6.01053.410

Das Metrosep C-Supp-2-Trennmaterial basiert auf einem Polystyrol-Divinylbenzol-Copolymer mit Carboxylgruppen. Es eignet sich zu Trennung und Bestimmung von ein- und zweiwertigen Kationen. Die Metrosep C Supp 2 - 100/4.0 Säule ist die kürzeste Trennsäule der Metrosep C-Supp-2-Familie. Sie eignet sich speziell zur Spurenanalytik der Standardkationen. Bestimmungsgrenzen unter dem µg/L-Bereich werden dank des extrem tiefen Basislinienrauschens nach sequentieller Suppression erreicht.

Lieferumfang 6.01053.410

Qt.	Order no.	Beschreibung
2 PCS	6.2744.060	Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.		



Optionales Zubehör

Order no.	Beschreibung
6.01053.500	Metrosep C Supp 2 Guard/4.0 The Metrosep C Supp 2 Guard/4.0 contains the C Supp 2 column material and is used to protect Metrosep C Supp 2 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The Metrosep C Supp 2 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.

