



Titrierkopf 6xNS14/ 3xNS9 (P&P)

6.01403.000

Titrierkopf mit 6 NS14 und 3 NS 9 Öffnungen für die Verwendung in Pick&Place-Modulen des OMNIS Sample Robot Pick&Place

Lieferumfang 6.01403.000

Qt.	Order no.	Beschreibung
2 PCS	6.1236.020	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm and O-ring.



5 PCS

6.1446.000

SGJ stopper / B-14/(15)



3 PCS

6.1446.010

Stopper

Used with Sample Changers



1 PCS

6.1446.160

Thermometer stopper, B-14/15, oblique insertion

To have the inserts of the 6.1458.010 titration head in a sloping position.



3 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders

