



Becheradapter zu OMNIS Probenbecher 120 mL

6.01404.030

Becheradapter zur Zentrierung des Probenbecher 120 mL im Schlitten der Pick&Place-Module

Optionales Zubehör

Order no.	Beschreibung
6.01400.200	Sample beakers, plastic (PP), 120 mL, 20 pieces
	Sample beakers, plastic (polypropylene, PP), 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.



6.01400.300 Sample beakers, clear glass, 120 mL, 20 pieces

Sample beakers, clear glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.



6.01400.303 Sample beakers, amber glass, 120 mL, 20 pieces

Sample beakers, amber glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces

