




Becheradapter zu OMNIS Probenbecher 75 mL

6.01404.040

Becheradapter zur Zentrierung des Probenbecher 75 mL im Schlitten der Pick&Place-Module

Optionales Zubehör

Order no.	Beschreibung	
6.01402.000	Sample beakers, clear glass, 75 mL, 30 pieces	
	Sample beakers, clear glass, 75 mL, suitable for 602041040 OMNIS sample rack, 30 pieces	

6.01402.003 Sample beakers, amber glass, 75 mL, 30 pieces

Sample beakers, amber glass, 75 mL, suitable for 602041040 OMNIS
sample rack, 30 pieces

