



OMNIS Probenrack, 9 x 250 mL, (PP)

6.02041.010

OMNIS Probenrack zu OMNIS Sample Robot Pick&Place passend für 9 Probenbecher. Folgende Probenbecher können verwendet werden: 6.01400.000, 6.01400.003, 6.01400.100.

Kunststoff: Polypropylen (PP)

Lieferumfang 6.02041.010

Qt.	Order no.	Beschreibung
1 PCS	6.02042.000	Holder for rack label (3x)

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Optionales Zubehör

Order no.	Beschreibung	
6.01400.000	Sample beakers, clear glass, 250 mL, 10 pieces Sample beakers, clear glass, 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.	
6.01400.003	Sample beakers, amber glass, 250 mL, 10 pieces Sample beakers, amber glass, 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.	
6.01400.100	Sample beakers, plastic (PP), 250 mL, 10 pieces Sample beakers, plastic (polypropylene, PP), 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.	