




OMNIS Probenrack, 16 x 120 mL, (PP)

6.02041.030

OMNIS Probenrack zu OMNIS Sample Robot Pick&Place passend für 16 Probenbecher. Folgende Probenbecher können verwendet werden: 6.01400.200, 6.01400.300, 6.01400.303.

Kunststoff: Polypropylen (PP)

Lieferumfang 6.02041.030

Qt.	Order no.	Beschreibung
1 PCS	6.02042.000	Holder for rack label (3x)

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Optionales Zubehör

Order no.	Beschreibung	
6.01400.200	Sample beakers, plastic (PP), 120 mL, 20 pieces Sample beakers, plastic (polypropylene, PP), 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.	
6.01400.300	Sample beakers, clear glass, 120 mL, 20 pieces Sample beakers, clear glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.	
6.01400.303	Sample beakers, amber glass, 120 mL, 20 pieces Sample beakers, amber glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces	
6.01404.030	Beaker adapter for OMNIS 120 mL sample beaker Beaker adapter for centering the 120 mL sample beaker in the slide of the Pick&Place modules	