




OMNIS Probenrack, 25 x 75 mL, (PP)

6.02041.040

OMNIS Probenrack zu OMNIS Sample Robot Pick&Place passend für 25 Probenbecher. Folgende Probenbecher können verwendet werden: 6.01402.000, 6.01402.003, 6.1459.400.

Kunststoff: Polypropylen (PP)

Lieferumfang 6.02041.040

Qt.	Order no.	Beschreibung
1 PCS	6.02042.000	Holder for rack label (3x)

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Optionales Zubehör

Order no.	Beschreibung	
6.01402.000	Sample beakers, clear glass, 75 mL, 30 pieces Sample beakers, clear glass, 75 mL, suitable for 602041040 OMNIS sample rack, 30 pieces	
6.01402.003	Sample beakers, amber glass, 75 mL, 30 pieces Sample beakers, amber glass, 75 mL, suitable for 602041040 OMNIS sample rack, 30 pieces	
6.01404.040	Beaker adapter for OMNIS 75 mL sample beaker Beaker adapter for centering the 75 mL sample beaker in the slide of the Pick&Place modules	
6.1459.400	Sample beaker, plastic (PE), 75 mL, 144 pieces Sample beaker, plastic (polyethylene, PE) 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, and 6.2041.800, suitable for 6.02041.040 OMNIS Software sample rack, 144 pieces.	