




Leitfähigkeitsmesszelle $c = 0.1 \text{ cm}^{-1}$ mit Pt1000 (Fixkabel)

6.0914.040

Leitfähigkeitsmesszelle aus Edelstahl mit Zellkonstante $c = 0.1 \text{ cm}^{-1}$, mit integriertem Pt1000-Temperaturfühler und Fixkabel (4 x B (4 mm)).

Dieser Sensor eignet sich für Messungen von tiefen Leitfähigkeiten (0 bis $300 \mu\text{S}/\text{cm}$) wie z.B. in deion. Wasser oder für Messungen nach USP <645>.

Der Sensor kann mit der Adapterbox (6.2103.160) am 856 Conductivity Module verwendet werden.

Lieferumfang 6.0914.040

Qt.	Order no.	Beschreibung
-----	-----------	--------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



Optionales Zubehör

Order no.	Beschreibung
-----------	--------------

6.1420.100 Flow-through measuring vessel
For 12 mm dia. electrodes (pH, metal, conductivity)



6.2324.010 Conductivity standard 100 μ S/cm 250 mL
Conductivity standard for calibration of conductivity measuring cells with cell constant = 0.1/cm.



6.2324.110 Conductivity standard 100 μ S/cm, 5 x 30 mL
Conductivity standard for calibration of conductivity measuring cells with cell constant = 0.1/cm.

