



Metrosep A Supp 5 - 250/2.0

6.1006.230

Die Metrosep A Supp 5 - 250/2.0 ist die Microbore- Hochleistungstrennsäule mit der sich selbst komplexe Trennprobleme einfach und reproduzierbar lösen lassen. Das Anwendungsspektrum für diese Säule geht weit über den Nachweis der Standardanionen hinaus. Die Metrosep A Supp 5 - 250/2.0 wird immer dort eingesetzt, wo höchste Trennleistung mit tiefsten Nachweisgrenzen und tiefem Eluentverbrauch kombiniert werden müssen.

Durch den geringen Eluentenfluss eignet sich diese Säule besonders für die IC-MS-Kopplung.

Optionales Zubehör

Order no.	Beschreibung
-----------	--------------

The Metrosep A Supp 5 Guard/2.0 reliably protects the 2 mm versions of the Metrosep A Supp 5 and Metrosep A Supp 7 IC anion columns against contamination from the sample or eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/2.0 highly recommended.



The Metrosep A Supp 5 S-Guard/2.0 reliably protects the 2 mm versions of the Metrosep A Supp 5 and Metrosep A Supp 7 IC anion columns against contamination from the sample or eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is connected directly onto the respective separation column via capillary connection. The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 S-Guard/2.0 highly recommended.

