



Metrosep A Supp 5 - 100/4.0

6.1006.510

Die Metrosep A Supp 5 - 100/4.0 erlaubt hocheffiziente und schnelle Trennungen. Diese Eigenschaft macht die Metrosep A Supp 5 - 100/4.0 zur Standardsäule für kurze Analysenzeiten und die Bestimmung spät eluierender Anionen (z.B. Perchlorat).

Optionales Zubehör

Order no.	Beschreibung
6.1006.500	Metrosep A Supp 5 Guard/4.0
	The Metrosep A Supp 5 Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contamination from the sample or the eluent. It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/4.0 highly recommended.



6.1006.540

Metrosep A Supp 5 S-Guard/4.0

The Metrosep A Supp 5 S-Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contaminations from the sample or eluent.



6.1020.520

Metrosep A Supp 10 Guard HC/4.0

The Metrosep A Supp 10 Guard HC/4.0 is the high-capacity version of the Metrosep A Supp 10 Guard/4.0.

