



Metrosep A Supp 5 - 150/4.0

6.1006.520

Die 150-mm-Variante der Metrosep A Supp 5 zeichnet sich durch sehr gute Trenneigenschaften aus. Hohe Bodenzahlen und ausgezeichnete Peaksymmetrien erleichtern das Arbeiten im unteren $\mu\text{g/L}$ -Bereich. Die Partikelgrösse von $5\text{ }\mu\text{m}$ trägt entscheidend zur Trennleistung dieser Säule bei. Die Metrosep A Supp 5 - 150/4.0 bietet die optimale Kombination von Selektivität und Kapazität, mit der sich selbst komplexe Trennaufgaben in kurzer Zeit sicher lösen lassen. Diese Eigenschaft macht die Metrosep A Supp 5 - 150/4.0 zu einer der besten universell einsetzbaren Standardsäulen der Ionenchromatographie.

Optionales Zubehör

Order no.	Beschreibung
-----------	--------------

The Metrosep A Supp 5 Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contamination from the sample or the eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/4.0 highly recommended.



The Metrosep A Supp 5 S-Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contaminations from the sample or eluent.

