



Metrosep A Supp 5 - 250/4.0

6.1006.530

Die Hochleistungstrennsäule aus dem Hause Metrohm mit extrem hohen Bodenzahlen für anspruchsvollste Trennaufgaben. Selbst komplexe Trennprobleme lassen sich mit der Metrosep A Supp 5 - 250/4.0 einfach und reproduzierbar lösen. Die hohe Kapazität der Säule erlaubt z. B. den Nachweis von 1 µg/L Bromat neben 150 mg/L Chlorid ohne Probenvorbereitung. Das Anwendungsspektrum für diese Säule geht weit über den Nachweis der Standardanionen hinaus. Die Metrosep A Supp 5 - 250/4.0 ist die Säule der Wahl wenn es darum geht, hohe Reinheitsstandards in der Halbleiterindustrie oder im Kesselspeisewasser von Kraftwerken sicher zu kontrollieren.

Optionales Zubehör

Order no.	Beschreibung
-----------	--------------

The Metrosep A Supp 5 Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contamination from the sample or the eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/4.0 highly recommended.



The Metrosep A Supp 5 S-Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contaminations from the sample or eluent.

