



Metrosep A Supp 5 - 50/4.0

6.1006.550

In weniger als 6 Minuten trennt die Metrosep A Supp 5 - 50/4.0 die sieben Standardanionen. Selbst Fluorid ist noch vom Injektionspeak getrennt und kann einwandfrei integriert werden. Die auf einem Polyvinylalkohol-Polymer basierende Säule zeichnet sich wie alle Säulen der A-Supp-5-Familie durch hohe Bodenzahlen und damit ausgezeichneter Trennleistung aus. Die Metrosep A Supp 5 - 50/4.0 ist die Säule der Wahl, wenn einfache Trennaufgaben in kurzer Zeit gelöst werden müssen – und zwar ohne Verzicht auf sehr niedrige Nachweisgrenzen.

Optionales Zubehör

Order no.	Beschreibung
-----------	--------------

The Metrosep A Supp 5 Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contamination from the sample or the eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/4.0 highly recommended.



The Metrosep A Supp 5 S-Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contaminations from the sample or eluent.

