



Metrosep C 6 - 100/4.0

6.1051.410

Die 100-mm-Version der Metrosep-C-6-Säule ist für die Bestimmung der Standardkationen beispielsweise in Trinkwasser vorgesehen. Es werden kurze Retentionszeiten erzielt, wobei Natrium und Ammonium trotzdem noch sehr gut getrennt werden. Die höhere Kapazität des C-6-Materials erlaubt höhere Probenvolumen.

Lieferumfang 6.1051.410

Qt.	Order no.	Beschreibung
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optionales Zubehör

Order no.	Beschreibung	
6.1051.500	Metrosep C 6 Guard/4.0	
	The Metrosep C 6 Guard/4.0 contains the C 6 column material and is used to protect against particles and contamination. This considerably lengthens the service life of the analytical separation column. The Metrosep C 6 Guard/4.0 works according to the "On Column Guard System" and is screwed directly onto the respective separation column with virtually no dead volume.	
6.1051.510	Metrosep C 6 S-Guard/4.0	
	The Metrosep C 6 S-Guard/4.0 contains the C 6 column material and is used to protect it against particles and contamination. This considerably lengthens the service life of the analytical separation column.	