



Elektrolytgefäß aus Kunststoff

6.1245.010

Elektrolytgefäß aus Kunststoff mit Keramikdiaphragma für Double-Junction-Referenzelektrode.

Optionales Zubehör

Order no.	Beschreibung
6.0728.110	Ag/AgCl reference electrode without electrolyte Inner reference system without electrolyte, with plug contact 5.5 mm, for use with 6.1245.010 electrolyte vessel made of plastic.



6.0728.120

Ag/AgCl reference electrode

Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.



6.0728.130

LL-Ag/AgCl reference electrode for CVS

Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.

