



KF-Adsorberrohr

6.1403.040

Zu KF-Ausrüstungen. Inklusive Deckel und O-Ring.

Optionales Zubehör

Order no.	Beschreibung
6.1414.030	KF titration vessel lid with 5 openings For volumetric KF titrations.



6.2730.010

Screw nipple

Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations.

