



Drehverschluss mit Septum, 100 Stk.

6.1448.067

Drehverschluss mit Septum, passend zu Probenglas 8 mL 6.2420.007, 100 Stück

Optionales Zubehör

Order no.	Beschreibung
6.1448.077	Septum, 100 pcs.

Replacement septum for 6.1448.067 screw cap, 100 pieces



6.2420.007

Sample glass with thread, 8 mL, 100 pieces

Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces


